

AIR-FREE CONDUCTIVE FABRIC BLANKETS & MATTRESSES



CLINICALLY SUPERIOR
PATIENT WARMING



IMPROVE PATIENT CARE WITH THESE BENEFITS:

SAFER

- Air-free warming.
- Prevents exposing patients to the risks associated with forced-air warming.
- Recommended by orthopedic surgeons during implant surgery for improving patient safety.

MORE EFFECTIVE

- Published studies show equal efficacy with forced-air warming.
- Customer experiences, however, often reveal superior efficacy.
- Capable of maintaining patient normothermia by warming from pre-op through PACU.

LESS EXPENSIVE

- No more costly disposables.
- Keep expenses flat and predictable with this easily cleanable, reusable warming system.
- Increase SCIP compliance, improve patient outcomes, and save money.



WHAT ARE USERS SAYING?

“Our unit has had a significant reduction in joint replacement infection since changing to an alternative patient warming system.”

-M.R. Reed, MD, Orthopedic Surgeon, Letters to the Editor, JBJS Feb 2011

AIR-FREE CONDUCTIVE FABRIC BLANKETS & MATTRESSES



THE NEXT WAVE
IN PATIENT WARMING

WHAT IS HOTDOG PATIENT WARMING?

Conductive fabric warming is state-of-the-art resistive, electric heating.

Here's how it works...

- HotDog controllers deliver low voltage electricity to the blankets and/or mattresses via thin, flexible cables, where it is converted into safe heat by the patented ThermAssure® conductive fabric.
- The lightweight, flexible heater fabric produces uniform and precisely controllable heat that is managed by a microprocessor in the controller and an advanced sensor system.
- The heater is sealed in an antimicrobial, non-porous shell--it's 100% cleanable.
- The result is safe and effective patient warming: an advancement in the mission to improve patient care.



CLINICALLY SUPERIOR PATIENT WARMING



VERSATILITY

- HotDog provides complete perioperative warming solutions, from pre-op through post-op.
- Whether for supine, lateral, or lithotomy positions or cardiac and burn units, HotDog is the warming solution for your most basic to your most challenging cases.



In addition to being:

- ☞ SAFER
- ✓ MORE EFFECTIVE
- 🐷 LESS EXPENSIVE

...HotDog is also:

- 🌿 GREEN
- 👤 ERGONOMIC

Learn more at hotdog-usa.com

About Augustine Temperature Management

Experience. Innovation. Dr. Scott Augustine revolutionized patient warming with the introduction of Bair Hugger®* forced-air warming 25 years ago. Since then, hundreds of millions of patients have received the benefits of normothermia during surgery.

Due to concerns with patient safety in ultra-clean surgeries and the spiraling costs that forced-air disposables put on the healthcare system, Dr. Augustine and his team of engineers sought to create an alternative warming solution to meet the needs of the healthcare community. The solution is HotDog air-free conductive fabric patient warming.

*Bair Hugger is a registered trademark of Arizant/3M.

Our Mission

We at Augustine Temperature Management are dedicated to improving patient care through the products we create. We aim to provide clinicians with a safer, more effective, and less expensive means of preventing the adverse effects of hypothermia during surgery. Patient safety is our passion.

www.hotdog-usa.com



6581 City West Parkway
Eden Prairie, MN 55344
+1 952-746-1720
1-888-439-2767

©2012 Augustine Temperature Management, LLC. HotDog is a registered trademark. All rights reserved. PN M120

